

SmartFvet™

# SmartFvet™

A Leading Veterinary Equipment Supplier

One Source For Vet Equipment ▶



# Smart Fvet

<https://smartfvet.com>



## CATALOGUE

1. Small Animal
2. Large Animal
3. Equine
4. Lab Diagnostics





Small Animal



Large Animal



Equine



Lab Diagnostics

# HV-Vet5 Anesthesia Workstation

10.2" TFT touch screen with adjusted angle, easier to observe and operate

Intelligent VS ventilation mode, minimum tidal volume 5ml



Multi-parameter Monitor shelf

Electrically driven and controlled, no drive gas needed

Intergraed Cycle Absorber with heating system Enhancing Comfort and Reduing Airway Damage

Vaporizer with compensation Compatible with Various Types of Vaporizers

Optional Medical-Grade Oxygen Concentrator, Bid farewell to bulky oxygen cylinders

Integrated Oxygen Concentrator Tray Design, Oxygen concentrator can be placed effortlessly, no need for fixation



# HV-X40

## Anesthesia Workstation

SmartFVet™

Anesthesia Workstation



PCV



SIMV



PSV(Optional)

### ❖ Features

- Suitable for pediatric and adults
- Tidal volume setting 50-1500mL
- Friendly user interface
- 8" LCD color screen with touch screen
- Knob and hard key Input

- **Sufficient modes of Ventilation:**VCV, PCV, SIMV-V, Manual
- **Gas supply:** O<sub>2</sub>, N<sub>2</sub>O, Air  
Optional: SIMV-P, CPAP/PSV, BACKUP, O<sub>2</sub>, ETCO<sub>2</sub>, AGSS
- 6 Tubes flowmeter more accurate and stable, auxiliary

#### Three gas supply with 6 tubes flowmeter

O<sub>2</sub>, N<sub>2</sub>O, Air, fresh supply, with 5 Gas Supply Pressure Gauges easy monitoring of central gas supplies and gas cylinders status.

#### Dual DIN rail sys

Adapt patient monitor support arm, BIS SYS or syringe pump frame on both sides (Optional).



#### Monitoring interface

Overall parameters of setting, measurement intuitive displays.



#### Anesthesia Gas Scavenging System

Independent active AGSS exhaust, avoid anesthesia air pollution in the OR, concerns clinicians healthy (Optional)



#### Antistatic casters

Diameter 125mm, two brakes of four casters.

#### Top shelf and light

Load weight 30Kg, Height: 1410 mm, Width: 550 Depth: 350,LED light illuminates all work surfaces.

#### Friendly interface

8.4" TFT touch screen, with rotary knob. Resolution 680 x 480



#### Anesthesia ventilator

Tidal volume: 50-1500mL. Easy to upgrade advance ventilation modes.

#### Dual Vaporizer

Select bar support 2 vaporizers position for Isoflurane, Sevoflurane, Halothane VP10 vaporizer compensates for variances in pressure, temperature and flow and needs no annual recalibration.



#### Yoke (REAR)

PISS standard, optional max 3 cylinders (Optional).



#### Drawer

Size: 200 x 392 x 398 (H x W x D) x2, Optional max 3.

#### Suction

Patient suction system with a regulator and reusable canister (Optional).




# HV-AM32


## Mobile Anesthesia Machine



Manometer



CO<sub>2</sub> Absorber



Breath Valve



Vaporizer


### Application

Small Animal  
Companion Anial



### Breathing Mode

Rebreathing Circuits  
Non-Rebreathing Circuits



# HV-1300


## Mobile Anesthesia Machine

SmartFVet™

Mobile Anesthesia Machine

### ❖ Features

- ★ Can be upgrade into CT & MR compatible model.
- ★ Adopts high accurate and reliable vaporizer.
- ★ Safety Pop-off valve design, enhance the safety and high efficiency.
- ★ Adopts module design, easy to maintain and repair.
- ★ The integrated circuit design greatly reduces breathing resistance and it is more suitable for animals.
- ★ Quick switch between rebreathing and non-rebreathing circuit, convenient to use.
- ★ Can be converted into a desktop model, suitable for more occasions.
- ★ Different attachment design, such as patient monitor stand, anesthesia gas scavengers, oxygen tank brackets to meet different customers' requirement



# HV-1300T

## Desktop Anesthesia Machine

SmartFVet™

### ❖ Features

- ★ Adopts high accurate and reliable performance vaporizer.
- ★ Safety Pop-off-Valve, enhance safety and high efficiency.
- ★ Adopts module design, easy maintain and repair.
- ★ Adopts integrated air path design, greatly reduce breathing resistance, more suitable for animals using.
- ★ Quickly switch between rebreathing and non-rebreathing circuit.
- ★ A variety of other equipment external hangers can be added to meet different needs of hospitals.






# HV-1300W

## Desktop Anesthesia Machine


SmartFVet™

Desktop Anesthesia Machine



# HV-1500

## Mobile Anesthesia Machine



### ❖ Features

- ★ Adopts high accurate and reliable performance vaporizer.
- ★ Safety Pop-off-Valve, enhance safety and high efficiency.
- ★ Adopts module design, easy maintain and repair.
- ★ Adopts integrated air path design, greatly reduce breathing resistance, more suitable for animals using.
- ★ Quickly switch between rebreathing and non-rebreathing circuit.
- ★ A variety of other equipment external hangers can be added to meet different needs of hospitals.



# HV-DAV80

## Anesthesia Ventilator

Bellow: special for small and big animal(0-300ml, 300-1500ml)

9 inch touch screen, set and monitor the breathing parameter

Silence key: press to mute for 120 seconds when alarming

Lock key: lock screen to prevent accidental touch

Shuttle knob: adjust different parameters



Breathing system interface  
The ventilator is connected to the CO<sub>2</sub> absorber

Power switch

Power cable interface

Gas source interface:  
Supply gas to ventilator

Driving gas outlet:  
supply gas to the bellows to make it work

Exhaust exhaust interface: connect the exhaust gas absorption tank to discharge exhaust gas

Driving gas inlet: connect the driving gas outlet to provide gas for bellows

Software upgrade usb interface

Built-in lithium battery, can be used for more than two hours

Flow sensor interface: monitor the flow parameter

Harmless exhaust gas outlet



# HV-JAY10 Oxygen Concentrator

## ❖ Features

- ★ Simple operation functions and easy-to-read LCD display.
- ★ Innovative cooling technology to ensure 24/7/365 operation.
- ★ Four flexible castors, easy to move or situate in any direction.
- ★ Humidifier bottle with a pressure relief value of 6PSI, safer and more comfortable.
- ★ Connected to HBOC, rapidly and effectively increase the patients' SPO<sub>2</sub>.
- ★ Connected to HFNC, provide oxygen with stable high concentration, more comfortable.
- ★ With medical grade oxygen concentration above 90%, better treatment effect for patients with pulmonary fibrosis.



# HV-GE5L Oxygen Concentrator

- ★ Set tired function via the display show the total work time
- ★ Set the sleep timer function, easy to use
- ★ Equipped with a barometric pressure of pressure relief valve, more security
- ★ Set the power off alarm, power outage in the course of their work reminds the user or caregiver
- ★ The design pressure, circulatory failure alarm function
- ★ The design of the compressor concentration alarm function
- ★ The design of a low oxygen concentration alarm function
- ★ The top of the design of the storage box can be placed on the attachment



# HV-PM12

## Multi-parameter Patient Monitor



- Silence
- Alarm Pause
- Freeze
- NIBP
- Record
- Main menu
- Encoder

12.1" Real Color TFT Screen

Built-in Lithium Battery

Standard Configuration

ECG, HR, RESP, SPO<sub>2</sub>, PR, NIBP, TEMP

Optional Features

ETCO<sub>2</sub>, Built-in Recorder

Six Major Functions



# HV-PM500K

## Multi-parameter Patient Monitor

Smart **F**Vet™

**[ Masimo oxygen saturation module from the United States ]**

Anti-motion, Anti-hypoperfusion

The "Gold Standard" in Oxygen Saturation Measurement



**SENSITIVITY**  
**97%** ↑

**SPECIFICITY**  
**95%** ↑

Alarm Performance real body motion and low per-fusion up to 97%

The low rate of false alarms body motion and low per-fusion 95%

 **SunTech Medical**<sup>®</sup>  
The Difference in Clinical Grade™

Built-in SunTech blood pressure module ensures more stable monitoring of animal blood pressure.



Mainstream end-tidal CO<sub>2</sub> module




Sidestream end-tidal CO<sub>2</sub> module



Wireless

OR



Wired

# HV-310

## Handheld Multi-gas Monitor

SmartFVet™

The HV310 Hand-held anesthesia machine monitor is a portable, user-friendly and compact device, measuring five anesthesia gas agents (ISO, SEV, ENF, DEF, HAL), N<sub>2</sub>O, ETCO<sub>2</sub>, FICO<sub>2</sub>, RR, SpO<sub>2</sub> and PR in ICU, OR and Emergency environment for Patients.

### ❖ Features

#### Display&resolution:

720P (1280 \* 720), 5.5 inch

#### Parameters Displayed:

ETCO<sub>2</sub>, FICO<sub>2</sub>, SPO<sub>2</sub>, RR, PR, AG

#### Indicators:

Battery status, Sound-light alarm,  
SpO<sub>2</sub> pulse pitch tone

**Interface:** USB 2.0

**Battery life:** > 8h

**IPower supply:**

90 -260V AC/5V DC, 50- 60 Hz, 20W

Built-in 8000mAh Li-ion battery



### Technical Specifications

Parameters	Range	Accuracy
ETCO <sub>2</sub>	0-15%	± (0.2% + 5% of reading)@0-10% ± (0.2% + 8% of reading)@10-15%
RR	0-150bpm	± 1bpm
SpO <sub>2</sub>	0-100%	± 2%@ (70- 100%)
PR	20-250bpm	± 2bpm
N <sub>2</sub> O	0-100	± (2% + 3% of reading)
Sevoflurane	0-10%	± (0.15% + 5% of reading)@0-5% ± (0.3% + 10% of reading)@5-10%
Isoflurane	0-8%	± (0.15% + 5% of reading)@0-5% ± (0.3% + 5% of reading)@5-8%
Enflurane	0-8%	± (0.15% + 5% of reading)@0-5% ± (0.3% + 5% of reading)@5-8%

# HV-F38

## Handheld Monitor



# HV-26

## Handheld Monitor



# HV-MT

## Handheld Monitor



# HV-300

## Portable EtCO<sub>2</sub> Capnography

Portable carbon dioxide monitor  
Measurements

EtCO<sub>2</sub> End-Tidal Carbon Dioxide

RR Respiration Rate





# HV-ECG3

## ECG Monitor

SmartFVet™

Complete data display

Full digital filtering, highresolution thermal printingPerfect ECG waveform

### ❖ ECG collection display

Input circuitFloating, Protection against defibrillator effect

Display Size: 7 "Display

Voltage tolerance:  $\pm 500\text{mV}$

Time constant: $\geq 3.2\text{s}$

Frequent response:  $0.05\text{Hz} \sim 165\text{Hz}$

Input Impedance:  $50\text{M} (10\text{Hz})$

Input Circuit Current:  $\leq 50\text{nA}$

Patient Auxiliary Current:  $< 10\text{A}$

( $220\text{V}-240\text{V}/50\text{Hz}$ )

CMRR: $> 60\text{dB}$

A/D ResolutionA/D: 12bits

Acquisition Mode:simultaneously 12 leads



# HV-MT02

## Blood Pressure Monitor



Basic model of blood pressure monitor PetPro

1. 4.3 "high resolution color LCD screen. 2. APP is available to measure and save data with mobile phone or tablet;
3. High accuracy in quiet and moving animals;
4. Standard BP cuffs in various sizes;
5. Multilingual: English, French, Spanish, Turkish, etc.;
6. Support 50 ID,each ID can store 200 data; Continuous monitoring can retain 1000 data;
- 7.Rechargeable lithium battery;
8. Bonus will be awarded to excellent dealer and top sales.

# HV-BLX9

## Dental X-RAY

### ❖ PRODUCT PARAMETER

#### Dental X-ray Machine

Power	16.8V/2600MA
Voltage	65KV
Current	0.2MA
Frequency	30KHz
Focus Point	0.3mm
Time Of Exposure	0.5-9s
Radiation Leakage	<20µm gy/h
Rate Power	60W



# HV-BLX6

## Dental X-RAY

### ❖ Features

- 1.The microprocessor-controlled timer brings high-tech X-ray imaging at your fingertips.
- 2.The control panel allows an accurate and fast exposure time selection while protecting the head during X-ray examinations.
- 3.Fitted with 60KV 8mA tube, the machine has been specially designed for intraoral dental diagnosis applications
- 4.Complies with local and international regulations on ionizing radiation protection for maximum safety
- 5.The soft positioner arm ensures naturally smooth and well-balanced movements,It's easy installed at the desired height and length.
6. It's ergonomics and design have been studied to meet the requirements of the most demanding practices in matters of floor space optimization.



# MD-480A Phaco Emulsifier

## ❖ Features



**Phaco Emulsification Frequency:**

40khz+20%

**Us Handpiece Top End Vibration Accuracy:**

Max.120m20%

**Stimulating Power Range:** 0-100%,

**Soft Keyboard Input Or Adjustable Stepping, Step:** 1%

**Pulse Frequency:** 1hz-99hz, Sot Keyboard Input Or Adjustable Stepping, Step:1hz;

**Efective Emulsification Time:** 1%-99% (E5ms)


**Soft Keyboard Input Or Adjustable Stepping: Step:** 1%;

**Display Accuracy: <+10%burst Time:** 5ms-100ms,

**Soft Keyboard Input Or Adjustable Stepping: Step:** 1ms:

**Display Accuracy: \$+10%max.vacuum:** 66.7kpa(500mmhg);

**Display Accuracy: S+15%max. Flow Rate:** 40 ml/min; Display Accuracy: S+20%





## HV-5S1 Ophthalmic Slit Lamp

### ❖ Features


- High Illumination White LED Light Source.
- Rechargeable Green Lithium Battery.
- Fully Charged Can Be Used for 5 or 6 Hours.
- Light-weight & Compact Design, Single Hand Operate.
- Link-able Eyepiece Equip with Four Spots.
- Can Be Digitized by Coopix Camera or Phone Eyepiece Adapters.



## HV-SL03 Ophthalmic Digital Slit Lamp

### ❖ Features

- High definition imaging
- Real-time focus system
- Super large touch screen
- Cloud service system



# HV-1000P

## Ophthalmic Fundus Camera

SmartFVet™

### ❖ Features

Dimensions: 30×9×19cm

Weight: 450g (typical)

Field of View Angle: 40±2°

Patient Refractive Compensation Range: -20D to +20D

Focus Method: Manual

Image Clarity (Image Pixels): 1920×1080Pixels

Infrared LED Wavelength Peak: 810nm±20nm

Flashlight Color Temperature: 4500K≤TC≤6700K

LCD Display: 3.5" FULL Color TFT LCD

Image Format: JPEG (Photograph)

File Transfer: Wi-Fi connection to iPad

File Storage: Micro SD Card, supports up to 16G

Power Types: Voltage ~220V, Frequency 50Hz, Lithium Battery 3.6V 3400mAh



# HV-2000P

## Ophthalmic Fundus Camera

### ❖ Features

★ @Non-mydratic and minimum pupil size 2.5mm

★ 10 Mega Pixels resolution

★ 5 targets fixation

★ auto capture

★ Andriod patient record managements system

★ 4.7" touch screen, with Wife, BT, and 4G transfer



# HV-D03

## Veterinary Endoscope

- ★ Support For Operation Video Playback. Network Transmission, Real-time Operation.
- ★ High Integrated: It Is Portable As Laptop Size Instead Of The Traditional Heavy Endoscope Tower, Convenient For Mobile Medical Treatment And Clinic Examination, Etc.
- ★ 17.3" Medical 4k Hd Monitor: 1.07 Billion Colors, Ultra-low Latency, Can Be Used As An Extended Display, Support Hdmi 4k Signal Input
- ★ LED Medical Light Source: Support Full Range Dimming, Illumination  $\geq 2000000$ lux
- ★ Camera Head: Ipx7 Waterproof, Supports Multiple Disinfection Methods Such As Ethylene Oxide And Alcohol, Compatible With Both Rigid Scopes And Flexible Scopes.



# HV-HD1

## Veterinary Endoscope

### ❖ Technical Parameters

**Resolution:** 3840\*2160P 60

**Touch screen:** 5inch LCD display

**Memory device:** High-resolution images and videos saved on USB storage

**Output:** DVI/HDMI/3G-SDI

**Language:** English/Chinese

**Product size:** 360mm\*300mm\*110mm

**LED light:** 100W Color index: Ra>92

**Light color temperature:** 5700K

**Bulb working life:** 20000h



# HV-1315P

## Veterinary Endoscope

### ❖ Features

- ★ Full HD (1080P) image processor for high-quality image visualization
- ★ Compatible with rigid endoscopes and flexible endoscopes
- ★ Multi levels LED light source and 3 levels air/water pump
- ★ 15.6" Full HD monitor with Medical Record Management System for better performance

SmartFVet™

# HV-E011

## Veterinary Endoscope

### ❖ Features

#### Multi-function Buttons:

1 . Zoom 2 . Capture/Record 3 . Freeze

**Divided structure :** Compatible with all products

**120° Left/Right Rotation :** Relieve fatigue of wrist rotation

**Video cable with 180° rotation :** Easy operation




# HV-E022

## Veterinary Endoscope

### ❖ Features

- ★ 1/2.8 SONY CMOS FULL HD Image Sensor
- ★ Built-in USB Video recording function;
- ★ Built-in LED Cold Light Source
- ★ The Integrated Design of LED Cold Light Source And Camer
- ★ IPX8 Waterproof Camera



## SV-T02 Animal Thermometer

- ★ Easy To Operate/Safe/Reliable/Hygienic
- ★ Automatic Shutdown/Beep Prompt
- ★ Quick Temperature Measurement
- ★ LCD Display/memory Function
- ★ Probe Lengthening/screen Enlargement
- ★ °C/°F Unit Convertible
- ★ Silent Measurement
- ★ Waterproof Type
- ★ Equipped With A Lanyard



## SV-T03 Animal Thermometer

### ❖ Features

- ★ Easy To Operate/Safe/Reliable/Hygienic
- ★ Automatic Shutdown/Beep Prompt
- ★ Quick Temperature Measurement
- ★ LCD Display/memory Function
- ★ Probe Lengthening/screen Enlargement
- ★ °C/°F Unit Convertible
- ★ Silent Measurement
- ★ Waterproof Type
- ★ Equipped With A Lanyard



## SV-T04 Animal Thermometer

### ❖ Features

- ★ Easy To Operate/Safe/Reliable/Hygienic
- ★ Automatic Shutdown/Beep Prompt
- ★ Quick Temperature Measurement
- ★ LCD Display/memory Function
- ★ Probe Lengthening/screen Enlargement
- ★ °C/°F Unit Convertible
- ★ Silent Measurement
- ★ Waterproof Type
- ★ Equipped With A Lanyard





# HV-402V

## Veterinary WIFI Otoscope

SmartFVet™

### ❖ Features

- ★ Exquisite and portable design.
- ★ fast and convenient way of focusing.
- ★ customized observation aperture (for ears).
- ★ one-touch for pictures and videos with 5x-80x.
- ★ Connected to Phones or PC by both WIFI and USB transmission.



# HV-S900 Veterinary Otoscope



## ❖ Features



**3" Color TFT IPS 640X480 LCD**  
True colors transmission



**4 super bright LEDs**  
for a brighter, whiter light

**VGA 640X480 pixels**  
Real-time observation



**LED dimming**  
Adjustable brightness

**Snapshot & Video**  
Micro SD up to 64GB



**Li-Ion rechargeable batteries**  
Come with a charger docking



**Connect to computers  
/ monitors**

# SV-C3LP

## Wireless Color Ultrasound



### ❖ How it works

The **SmartFVet**™ Wireless Probe Is A Mini Ultrasound Scanner Without A Screen. We Minimize The Components Of A Traditional Ultrasound Into A Small Circuit Board Built Into The Probe, And Showing Image In Smart Phone/tablet Through Wifi Transferring. rotate Built-in Screen Probe Is A Wireless Ultrasound Probe Add Screen, Image Can Both Show In Screen And Tablet. Image Transferring Through Internal Wifi From Probe, No Need External Wifi Signal.



# SVProbe

## Wireless Color Ultrasound

### ❖ Features

- ★ Advanced digital imaging technology, clear image
- ★ High cost-effective
- ★ Wireless connectivity, wireless charging, easy to operate
- ★ Small and light, easy to carry
- ★ Whole probe is waterproof, convenient for cleaning.



# SUprobe

## B Wireless Ultrasound



- Display:** iPad Series
- Display mode:** B mode
- Frequency:** 5MHz
- Gray Scale:** 256
- Scanning mode:** Linear Array
- Size:** 145mm\*45mm\*40mm
- Scanning Depth:** 40-120mm
- Battery last:** 4 hours
- Tranducer Length:** 45mm
- Weight:** <300g

### How it works?

The Sonostar wireless probe is a mini ultrasound scanner without a screen. We minimize the components of a traditional ultrasound into a small circuit board built into the probe, and showing image in smart phone/tablet through Wifi transferring. image can both show in screen and tablet. Image transferring through internal wifi from probe, no need external Wifi signal. Only one probe, connect to tablet by wifi, convenient for use. And with special function for backfat measure function.

# SV5 Palm Ultrasound

## ❖ Feature

- ★ 7 inch TFT screen
- ★ Pretty and very light
- ★ Full digital imaging technology, crystal-clear Image
- ★ Can work with backfat probe
- ★ Low power consumption
- ★ 2G Large volume storage image
- ★ Built-in battery
- ★ Solid shell and structure
- ★ Waterproof(Optional)
- ★ High Cost-effective



# HV-20A Palm Ultrasound

## ❖ Feature

- ★ Handheld , light weight
- ★ Should & waist belts allow free hands working
- ★ Rubber protector around the scanner against dirt and impact from the outside
- ★ Splash water
- ★ Both battery and AC powered




# HV-C5Pro

## Portable Ultrasound

SmartFVet™

### ❖ Specification

- ★ **Screen** : 15inch LCD
- ★ **Image mode** : B, B/M, Color, PW, PDI
- ★ **Channel of Ultrasound**: 64
- ★ **Gain** : 30db-105db
- ★ **Denoise** : 0-1-2-3-4
- ★ **Dynamic Range** : 40-50-60-70-80-90-100-110
- ★ **Probe type** : Convex, Linear, Phase Array, Micro-Convex, Transvaginal, Transrectal
- ★ **Probe Connection Port**: 2
- ★ **Frequency, Scanning Depth, Angle/Width**:
  - Convex : 3.5MHz/5MHz,90/160/220/305mm,60°
  - Linear : 7.5MHz/10MHz, 20/40/60/100mm,40mm
  - High Frequency Linear : 10MHz/14MHz,20/30/40/55mm,25mm
  - Phase Array : 6.5MHz/8MHz,40-100mm,80°
  - Micro-Convex : 3.5MHz/5MHz,90/130/160/200mm,88°
  - Transvaginal : 6.5MHz/8MHz,40/60/80/100mm.149°
  - Biplane Transrectal : Convex 6.5MHz/8MHz40/60/80/100mm, 149°; Linear 7.5MHz/10MHz,20/40/60/100mm,40mm
- ★ **Image Adjust** : BGain, TGC, DYN, Focus, Depth,Harmonic,Denoise, Color Gain, Steer, PRF
- ★ **Cineplay**: auto and manual, frames can set as100/200/500/1000
- ★ **Puncture Auxiliary Function**: In-Plane puncture guide line, Out-Plane puncture guide line,automatic blood vessel measurement.
- ★ **Measure**: Length, Area, Angle, Heart Rate, Obstetrics- Image Frame Rate: 18f/s
- ★ **Operating System** : Windows



## SV-P10

### Palm Ultrasound

- ★ 15 inches LCD display zone for large image size, with high quality.
- ★ Probe Interface: Dual probe interface, with automatic identification function
- ★ ARM platform, fast computing, HD imaging
- ★ Abundant probe sets, suitable for different parts.
- ★ Battery Capacity: Independent battery compartment /removable lithium battery /5200 mAh
- ★ Four encoders, knob design, quick and convenient operation
- ★ Support multi-language switching
- ★ Movie playback: 256 frames, can be played back frame by frame and pause playback, pause playback can be used for image preservation and measurement
- ★ Measurement: Professional measurement software;



# HV-MRI1.5T

## MRI

SmartFVet™

### ❖ Feature

- ★ Superior Images
- ★ Low Total Cost Of Ownership
- ★ Easy To Use
- ★ Head Coil
- ★ V-shaped Patient Bed And Coil
- ★ Flat Patient Bed And Flex Coils
- ★ Integrated Rf Shield



### The operator room system includes

- ★ Monitor, keyboard, mouse and USB connection
- ★ Workstation software

### ❖ Superconducting magnet

Actively shielded superconducting magnet system with integrated closed cycle refrigerator system:

**Field strength:** 1.5 Tesla

**Polarity:** North pole at patient end

**Homogeneity (shimmed):** 50cm dsv < +/-10ppm pk to pk  
45cm dsv < +/-5ppm pk to pk  
30cm dsv < +/-1ppm pk to pk  
Patient bore diameter: 60cm

**Magnet length:** 1.73m

**Stray field:** < 3.8m (Axial), < 2.5m (Radial)

### The equipment room subsystem consists of 4 major cabinets:

- ★ RF cabinet, including:
  - ★ 8 channel spectrometer, RF amplifier, Gradient controller, Equipment room computer, Magnet monitoring unit
- ★ Gradient cabinet
- ★ Power distribution unit
- ★ Water cooler cabinet

# HV-2012D-1

Ceiling Surgery Lighting



# HV-2036B-2

Wallmount Lighting



# HV-2012B-1

Examination Light



# HV-2036L-3 Examination Light

Veterinary Surgical Lighting

SmartFvet™



Lighting



# HV-2012L-1 Examination Light



# HV-IP01 Infusion Pump

## ❖ Features

- 4.3 inch color touch screen for easy operation and clear vision
- Compact and lightweight design for easy carrying and mounting
- User friendly interface with audible and visual alarm system
- Built-in lithium battery, around 5 hours working time after fully charged
- Multiple infusion modes to satisfy different infusion requirements
- Advanced pressure detection technology



# HV-IP02 Infusion Pump

# HV-MI20 Infusion Pump

## ❖ Features



# HV-SP01

## Syringe Pump

### ❖ Features

- 4.3 inch color touch screen for clear vision
- Simple rate mode with one key start , easy to operate
- Compatible with various brands of syringe
- Various infusion modes to meet different clinical requirements
- User friendly interface with audible and visual alarm system
- Rechargeable lithium battery can work continuously for around 5 hours

SmartFVet™



Pumps, Poles And Warmers



# HV-SP02

## Syringe Pump

# HV-MS51

## Syringe Pump



# HV-50 Syringe Pump

SmartFVet™

## ❖ Features

3.5" iTouch screen operation

**Syringe size :** 5ml,10ml,20ml, 30ml, 50ml,60ml.

Optional built-in WiFi module.

Can connected to HIS/CIS via WLAN or LAN.

Syringe slider protection

Internal battery up to 20 hours continuous powersupply.



# HV-6010 Infusion Pump

## ❖ Features

**Accurate:** A compact peristaltic pump structure is adopted to ensure the stability and accuracy of infusion.

**Easy to Use:** Outstanding touchscreen for fast responses to all settings and operation.

**Multiple infusion modes:** Meet different section and drug infusion requirements. Convenient to select.

### **Automatic pump door mechanism:**

5. Press lightly and the door opens. Patently designed automatic door guarantees infusion safety under complete closed pump door.

### **Reliable:**

In compliance with official guidance, SYS-6010 VET infusion pump integrates strict risk management control.



# HV-ICU05

## Icu incubator

### ❖ Features

1. Accurate Temperature Control--Being Equal Everywhere
2. Accurate Humidity Control—Healthy and Balanced
3. Negative-ion Generation—Truly Auxiliary Medical Effect
4. Sterilization Function ---- Free from Cross-contamination
5. Medical Atomization Treatment Function---- Making the Treatment More Convenient
6. Carbon Dioxide Concentration Monitoring Systems----A Defense Line to Guard Life
7. ICU Illumination Function---- Creating Comfortable Therapeutic Environment
8. Setting of Security Mechanism---- No Worries in Use
9. Human Oriented Structural Design ---- Reasonable Operation, Clean and Convenient



# HV-WA01

## Air Warming

### ❖ Material

Two sensors are located at the out let from the heat unit also at the end of the hose, detecting and compensating for heat loss from the hose.

### ❖ Features

- \* The operating noise level is very low. (53DBA).
- \* Accurate LCD display, real time temperature of blanket monitoring.
- \* Swivelling easy to fit outlet port, makes attaching the hose effortless.

# HV-210V

## Veterinary Warming System

SmartFVet™

Veterinary Warming System utilizes the convection of air to transmit the hot air generated by the main unit to the through the pipeline, so that the heat energy is evenly distributed and transmitted to the blanket to achieve heat preservation or heating effect. For the prevention and treatment of hypothermia in small animals, it can accelerate the recovery of small animals and shorten the hospitalization time of small animals.

- ★ Multiple Wind Speed.
- ★ Multiple Temperature Control Protection.
- ★ History Record.
- ★ Air Purification Function.





# HV-V5

## Veterinary Full Auto Dry chemistry Analyzer

**Features:** Up to 17 reportable parameters in a test, 32 parameters configured into 9 profiles.

**Data Input Display:** 8.0 inch capacitive touchscreen

**Print Report:** Built in thermal printer, External Printer via MNCHIP Medical Data Management Platform.

**Mode of Operation:** Continuous

**Ambient Temperature:** 10~30°C (50~86 °F), indoor use.

**Power Requirements:** 120 VA

**Main Supply Voltage:** 100-240 volts AC, 50-60 Hz.

**Humidity:** 40%~85%

**Reaction Temperature:** 37°C ± 0.2°C

**Reagent disc:** disposable, self-contained reagent in dry form




# HV-120VP


## Chemistry Analyzer



Small and portable




Easy to use



Little sample volume



One analyzer for multiple uses



Consumable advantages



Minimal maintenance



# SV-TG10 Centrifuge

SmartFVet™

## ❖ Features

- \* The CombiSpin Centrifuge is ergonomically designed for centrifuging relativelySilent, constant and stable, touching panel.
- \* Small volumes of biological fluids at speeds of up to 10,000 rpm.  
It is particularly suited to Veterinary Practice laboratories.
- \* The CombiSpin has been designed for fast, simple, intuitive and safe operation.
- \* The single rotor can hold both microhaematocrit capillaries and larger plastic tubes up to 2ml in volume therefore avoiding the need to change rotors.
- \* Three pre-set optimised speed/time programs are available for serum/plasma, microhaematocrits, and urine sediments and accessed by one click of a key.
- \* However, both speed and time can each be simply changed incrementally in case different settings are required for 'special' samples.
- \* Plastic tubes of up to 2ml in volume are placed in angled tube holders in the rotor and microhaematocrit glass capillaries are fitted horizontally in grooves on top of the rotor.
- \* A screw on rotor lid is provided for use when running microhaematocrits (it is not needed when running plastic tubes on their own).





# HV-T101

## Diode Laser

### ❖ Features

### You can choose:

- ★ Single dental treatment configurations.
- ★ Single physiotherapy configuration,
- ★ Both Together.




# HV-350D Electrosurgical Unit

SmartFVet™

## ❖ Features

- Remote fault diagnosis.
- Human voice alarm function.
- Full-function automatic power compensation.
- System has return electrode monitoring (REM) system.
- Microprocessor control and digital display in LED panel.
- Used for Monopolar and Bipolar surgery in many surgical fields.
- Support bipolar closure and separation to meet the requirements of surgical bipolar treatment without sticking the negative plate.
- System has a universal endoscope interface, which is used with gastroscope, enteroscope, bronchoscope, aparoscope and other endoscopes.
- This machine has a power-off protection circuit, which can memorize the output set value of each function in real time.
- Monopolar output can be controlled either by hand or foot, and monopolar/bipolar automatic output.
- Patient circuit bipolar electrode plate contact quality detection system.
- Support in underwater cut operation.
- Passed EMC test.



# HV-ZSG1 Ultrasonic Scalpel System

## ❖ Stable and durable performance

- ★ 55,000 calculations per second are performed, and the maximum conversion ratio of electric energy to vibration efficiency is always maintained, so as to minimize the extra heat and ensure the stability of the system working for a long time
- ★ The cutter head is continuously excited for more than eight hours without cutting





**SmartFvet**<sup>TM</sup>

One source for vet equipment

Equipment is curated for quality

Experience the difference

Competitive OEM/ODM factory

Full customization

Room A206, 2 / F, Building 1, Zhongdian Jinyang Science and Technology Park, No. 1, South Disheng Street, Economic and Technological Development Zone, Beijing, China

Web Site:<https://smartfvet.com>

Email:[info@smartfvet.com](mailto:info@smartfvet.com)

Phone:+86-199-4372-0031

WhatsApp:+86-199-4372-0031