



SmartFvetTM

Professional Support for Veterinarians
Equipment is curated for quality
Experience the difference
Competitive OEM/ODM factory
Full customization

Room A206, 2 / F, Building 1, Zhongdian Jinyang Science and Technology Park, No. 1, South Disheng Street, Economic and Technological Development Zone, Beijing, China

Web Site:<https://smartfvet.com>

Email:info@smartfvet.com

Phone:+86-199-4372-0031


WhatsApp:+86-199-4372-0031

SmartFvetTM

SmartFvetTM

A Leading Veterinary Equipment Supplier

Professional Support for veterinarians ▶



HV-HD1

Veterinary 4K Endoscope Camera System

SmartFVet™

HV-HD1 is a 4K endoscope camera that produces sharp and clear image, ideal for applications in various specialties. With the latest image technology, the camera helps surgeons in diagnosis and surgery with fine details and true color that will not be seen using a HD endoscope camera. The system has also been redesigned to improve operability for surgeons.

❖ Features

- ★ Producing native 4K image that is four times the resolution of HD
- ★ ENDO mode preset for specialization
- ★ A wider color gamut that ensures better observation for surgeons.
- ★ Compact and lightweight
- ★ Water-resistant
- ★ Programmable buttons

❖ Technical Parameters

Resolution: 3840*2160P 60

Touch screen: 5 inch LCD display

Memory device: High-resolution images and videos saved on USB storage

Output: DVI/HDMI/3G-SDI

Language: English/Chinese

Product size: 360mm*300mm*110mm

LED light: 100W

Color index: Ra>92

Light color temperature: 5700K

Bulb working life: 20000h



HV-D03

Veterinary Endoscope

SmartFVet™

❖ Features

- ★ Support For Operation Video Playback. Network Transmission, Real-time Operation.
- ★ High Integrated: It Is Portable As Laptop Size Instead Of The Traditional Heavey Endoscope Tower.
Convenient For Mobile Medical Treatment And Clinic Examination, Etc.
- ★ 17.3" Medical 4k Hd Monitor: 1.07 Billion Colors, Ultra-low Latency, Can Be Used As An Extended Display, Support Hdmi 4k Signal Input
- ★ LED Medical Light Source: Support Full Range Dimming, Illumination ≥ 2000000 lux
- ★ Camera Head: Ipx7 Waterproof, Supports Multiple Disinfection Methods Such As Ethylene Oxide And Alcohol, Compatible With Both Rigid Scopes And Flexible Scopes.



HV-E011

Portable Veterinary Endoscope System

SmartFVet™

❖ Features

★ Multi-function Buttons

①Zoom ②Capture/Record ③Freeze

★ 120° Left/Right Rotation

Relieve fatigue of wrist rotation

★ Detachable Structure

Compatible with all products

★ Video Cable With 180° Rotation

Easy operation

processor with 10 inch touch-screen.



HV-E022

Veterinary HD Endoscope

SmartFVet™

❖ Features

★ 1/2.8 SONY CMOS FULL HD Image Sensor

★ Built-in USB Video recording function


★ Built-in LED Cold Light Source

★ The Integrated Design of LED Cold Light Source And Camer

★ IPX8 Waterproof Camera



1920*1080P
SONY CMOS HD Camera



HV-1315P

Veterinary Endoscope Camera System

SmartFVet™

❖ Features

- ★ Full HD (1080P) image processor for high-quality image visualization
- ★ Compatible with rigid endoscopes and flexible endoscopes
- ★ Multi levels LED light source and 3 levels air/water pump
- ★ 15.6" Full HD monitor with Medical Record Management System for better performance



HV-602K-2

Veterinary Endoscope

SmartFVet™

❖ Features

- ★ Software Placement: Private computers Turn into Workstations
- ★ Micro Light Source: No Fiber Docking
- ★ **Concentration is the essence**

额定电压	灯泡功率	显色指数	靶面尺寸	有效像素	传感器类型	最大分辨率	焦距
Rated voltage	Bulb power	CRI	Target size	Effective pixels	Sensortype	Maximum resolution	Focal length
AC90~240V 50Hz~60Hz	LED 3W	≥90%	1/1.8"	6,300,000	CMOS	3088X2064	14~35mm

