



SmartF VetTM

Professional Support for Veterinarians
Equipment is curated for quality
Experience the difference
Competitive OEM/ODM factory
Full customization

Room A206, 2 / F, Building 1, Zhongdian Jinyang Science and Technology Park, No. 1, South Disheng Street, Economic and Technological Development Zone, Beijing, China

Web Site:<https://smartfvet.com>

Email:info@smartfvet.com

Phone:+86-199-4372-0031

WhatsAPP:+86-199-4372-0031

SmartF VetTM

SmartF VetTM

A Leading Veterinary Equipment Supplier

Professional Support for veterinarians ▶



smartFvet Monitors CATALOGUE


Multi-parameter Monitor	2
Multi-gas Monitor	4
Handheld Monitor	5
Portable EtCO2 Capnography	7
ECG Monitor	8
Blood Pressure Monitor	9



HV-PM12 Multi-parameter Patient Monitor



Multi-parameter Monitor



12.1" Real Color TFT Screen

Built-in Lithium Battery

Standard Configuration

ECG, HR, RESP, SPO₂, PR, NIBP, TEMP

Optional Features

ETCO₂, Built-in Recorder

Six Major Functions



HV-PM500K

Multi-parameter Patient Monitor



[Masimo oxygen saturation module from the United States]

Anti-motion, Anti-hypoperfusion

The "Gold Standard" in Oxygen Saturation Measurement



SENSITIVITY
97% ↑

SPECIFICITY
95% ↑

Alarm Performance real body motion and low per-fusion up to 97%

The low rate of false alarms body motion and low per-fusion 95%



Built-in SunTech blood pressure module ensures more stable monitoring of animal blood pressure.




Mainstream end-tidal CO₂ module

Sidestream end-tidal CO₂ module



Wireless

OR



Wired

HV-310

Handheld Multi-gas Monitor



The HV310 Hand-held anesthesia machine monitor is a portable, user-friendly and compact device, measuring five anesthesia gas agents (ISO, SEV, ENF, DEF, HAL), N₂O, ETCO₂, FICO₂, RR, SpO₂ and PR in ICU, OR and Emergency environment for Patients.

◆ Features

Display&resolution:

720P (1280 * 720), 5.5 inch

Parameters Displayed:

ETCO₂, FICO₂, SPO₂, RR, PR, AG

Indicators:

Battery status, Sound-light alarm,

SpO₂ pulse pitch tone

Interface: USB 2.0

Battery life: > 8h

IPower supply:

90 -260V AC/5V DC, 50- 60 Hz, 20W

Built-in 8000mAh Li-ion battery



Technical Specifications

Parameters	Range	Accuracy
ETCO ₂	0-15%	± (0.2% + 5% of reading)@0-10% ± (0.2% + 8% of reading)@10-15%
RR	0-1 50bpm	± 1bpm
SpO ₂	0-1 00%	± 2%@ (70- 100%)
PR	20-250bpm	± 2bpm
N ₂ O	0-100	± (2% + 3% of reading)
Sevoflurane	0-10%	± (0.15% + 5% of reading)@0-5% ± (0.3% + 10% of reading)@5-10%
Isoflurane	0-8%	± (0.15% + 5% of reading)@0-5% ± (0.3% + 5% of reading)@5-8%
Enflurane	0-8%	± (0.15% + 5% of reading)@0-5% ± (0.3% + 5% of reading)@5-8%

HV-F38

Handheld Monitor

SmartFVet™



HV-26

Handheld Monitor



SmartFVet™

HV-MT

Handheld Pulse Oximeter Monitor

◆ Features

- ★ The smallest and exquisite handheld vital sign monitor in the worldwide.
- ★ Accurate measurement in 0.05% Low perfusion Index, the world leading production in the movement environment.
- ★ Provide the parameter of PI (Perfusion Index 0.05%-20%).
- ★ Accurate measurement in all kind of extreme environment measurement such as neonatal and shock patients.
- ★ Multiple vital signs measurement parameters for collocation optional:
 - NIBP, SpO₂, TEMP
 - NIBP, TEMP
 - SpO₂, TEMP
 - SpO₂, CO₂
- ★ two measurements mode: Spot, Monitoring.
- ★ Data storage capabilities, support for historical data view and trend chart back to see.
 - In monitoring mode, can continuously record 48 hours of measurement data.
 - In spot mode, 100 users ID, each user can store 200 sets data.
- ★ Support audio and visual alarm,
- ★ support to set alarm upper and lower limit values.
- ★ Lithium ion battery power supply. Optional Charging dock.



HV-300

Portable EtCO₂ Capnography

Portable carbon dioxide monitor

Measurements

EtCO₂ End-Tidal Carbon Dioxide

RR Respiration Rate



HV-ECG3

ECG Monitor



Complete data display

Full digital filtering, high resolution thermal printing Perfect ECG waveform

❖ ECG collection display

Input circuit Floating, Protection against defibrillator effect

Display Size: 7" Display

Voltage tolerance: ±500mV

Time constant: ≥3.2s

Frequent response: 0.05Hz ~ 165Hz

Input Impedance: 50M (10Hz)

Input Circuit Current: ≤50nA



HV-UN8003V

ECG Monitor

❖ Features

- ★ Compact in size and ergonomic design, very easy to operate, portable.
- ★ 80mm, three channel format recording with interpretation;
- ★ Recording format: 1ch+V, 3ch, 3ch+V.
- ★ Three working modes: Auto, Manual and Arrhythmia Analysis.
- ★ High accuracy digital filter; Automatic baseline adjustment.
- ★ 7 leads simultaneously collection and display.
- ★ 480x272 graphic 3.5 inch color LCD displays detailed ECG information.
- ★ 100 cases of ECG can be stored and displayed.
- ★ Adapt to 110-230V, 50/60Hz power supply.
- ★ Rechargeable Ni-MH battery supports continuous work over 3 hours.
- ★ With Rs232 and USB communication interface, optional for connection to Laser.
- ★ PC ECG Management software is optional.



HV-MT02

Blood Pressure Monitor



❖ Features

- Basic model of blood pressure monitor PetPro
1. 4.3 "high resolution color LCD screen. 2. APP is available to measure and save data with mobile phone or tablet;
 3. High accuracy in quiet and moving animals;
 4. Standard BP cuffs in various sizes;
 5. Multilingual: English, French, Spanish, Turkish, etc.;
 6. Support 50 ID, each ID can store 200 data; Continuous monitoring can retain 1000 data;
 7. Rechargeable lithium battery;
 8. Bonus will be awarded to excellent dealer and top sales.



Three color options

